

LAB REPORT

Batch: 2605

Date:

	ACTIONS	RECORDED TEMP.	NOTES
DAY 1 . 0h - 24h			
DAY 2 . 24h - 48h	Turn at: 48 hrs		
DAY 3 . 48h - 72h			
DAY 4 . 72h - 96h	Turn at: 96 hrs	42	
DAY 5 . 96h - 120h		45.5	
DAY 6 . 120 - 144h	Move to dryer at 144 hrs	40.3	
DAY 7 . 144 - 168h			

NOTES: