

# LAB REPORT

Batch: 2604

Date:

	ACTIONS	RECORDED TEMP.	NOTES
<b>DAY 1</b> . 0h - 24h			
<b>DAY 2</b> . 24h - 48h	Turn at: 48 hrs		
<b>DAY 3</b> . 48h - 72h			
<b>DAY 4</b> . 72h - 96h	Turn at: 96 hrs	41	
<b>DAY 5</b> . 96h - 120h		48	
<b>DAY 6</b> . 120 - 144h	Move to dryer at 144 hrs	44	
<b>DAY 7</b> . 144 - 168h			

**NOTES:**